

# Anti-Proteasome subunit alpha type 6 Rabbit antibody



Catalog Number: 167768

## Product name

Anti-Proteasome subunit alpha type 6 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Proteasome subunit alpha type 6 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Proteasome subunit alpha type 6

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000-1/5000

## Validations



Western blot detection of Proteasome subunit alpha type 6 in A549,HL-60 cell lysates using Proteasome subunit alpha type 6 Rabbit mAb(1:1000 diluted).Predicted band size:27kDa.Observed band size:27kDa.