

# KLF13 antibody



Catalog Number: 112078

## Product name

KLF13 antibody

## Specificity

Human; other species not tested.

## Antibody description

KLF13 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HepG2 cells, K-562 cells. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 38-45 kDa

## Preparation

This antibody was obtained by immunization of KLF13 recombinant protein (Accession Number: NM\_015995). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:112078(KLF13 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human brain using Catalog No:112078(KLF13 antibody) at dilution of 1:50 (under 10x lens)