

Anti-MOS antibody

Catalog Number: 176078

Product name

Anti-MOS antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to MOS

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mos. aa range:61-110

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western Blot analysis of various cells using Mos Polyclonal Antibody