

# C10orf62 Polyclonal Antibody



Catalog Number: 163478

## Product name

C10orf62 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to C10orf62

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human C10orf62 (NP\_001009997.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - C10orf62 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat testis using C10orf62 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - C10orf62 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse testis using C10orf62 antibody at dilution of 1:100 (40x lens).