

PACSIN1 antibody



Catalog Number: 113478

Product name

PACSIN1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PACSIN1 Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells. Positive WB detected in mouse cerebellum tissue, mouse brain tissue. Observed molecular weight by Western-blot: 48-51kd

Preparation

This antibody was obtained by immunization of PACSIN1 recombinant protein (Accession Number: NM_020804). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:113478(PACSIN1 antibody) at dilution of 1:100



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:113478(PACSIN1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)