

# RASGRP2 antibody



Catalog Number: 114478

## Product name

RASGRP2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RASGRP2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse liver tissue, mouse lung tissue. Observed molecular weight by Western-blot: 48kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_153819). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114478(RASGRP2 antibody) at dilution of 1:300