

SNX25 antibody

Catalog Number: 115478

Product name

SNX25 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SNX25 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse brain tissue, mouse thymus tissue, Raji cells, rat brain tissue. Observed molecular weight by Western-blot: 98-105 kDa, 65 kDa

Preparation

This antibody was obtained by immunization of SNX25 recombinant protein (Accession Number: BC029868). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115478(SNX25 antibody) at dilution of 1:500