

TTI1 antibody

Catalog Number: 116478

Product name

TTI1 antibody

Specificity

Human; other species not tested.

Antibody description

TTI1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of TTI1 recombinant protein (Accession Number: XM_047440606). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116478(KIAA0406 Antibody) at dilution of



1:1000

IP Result of anti-KIAA0406 (IP:Catalog No:116478, 4ug; Detection:Catalog No:116478 1:1000) with HeLa cells lysate 1600ug.