

PMS2 antibody



Catalog Number: 107478

Product name

PMS2 antibody

Immunogen

[Human PMS2 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

PMS2 Mouse Monoclonal antibody. Positive WB detected in A431 cells, HeLa cells. Positive IHC detected in human cervical cancer tissue. Observed molecular weight by Western-blot: 100kd

Preparation

This antibody was obtained by immunization of PMS2 recombinant protein (Accession Number: NM_000535). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:107478(PMS2 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:107478(PMS2 Antibody) at dilution of 1:50 (under 40x lens)