

# GATA1 antibody



Catalog Number: 110878

## Product name

GATA1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GATA1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HL-60 cells, NIH/3T3 cells, Raji cells. Positive IHC detected in mouse spleen tissue. Observed molecular weight by Western-blot: 55kd

## Preparation

This antibody was obtained by immunization of GATA1 recombinant protein (Accession Number: BC009797). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110878(GATA1 antibody) at dilution of 1:1200



Immunohistochemistry of paraffin-embedded mouse spleen slide using Catalog No:110878(GATA1 Antibody) at dilution of 1:50