

Jagged 2 antibody



Catalog Number: 111878

Product name

Polyclonal

Jagged 2 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

Jagged 2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 150 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002226). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111878(JAG2 Antibody) at dilution of 1:500

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality