

XRN2 antibody

Catalog Number: 116878

Product name

XRN2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

XRN2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, COLO 320 cells, HeLa cells. Observed molecular weight by Western-blot: 104kd

Preparation

This antibody was obtained by immunization of XRN2 recombinant protein (Accession Number: BC006417). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116878(XRN2 antibody) at dilution of 1:200