

DBNL antibody

Catalog Number: 109878

Product name

DBNL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DBNL Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, human heart tissue, Jurkat cells, K-562 cells, mouse heart tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human breast cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of DBNL recombinant protein (Accession Number: NM_001014436). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109878(DBNL antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109878(DBNL antibody) at dilution of 1:100 (under 10x lens)


DBNL antibody

Catalog Number: 109878



Immunofluorescent analysis of HepG2 cells, using DBNL antibody Catalog No:109878 at 1:50 dilution and FITC-labeled donkey anti-rabbit

IgG(green).



IP Result of anti-DBNL (IP:Catalog No:109878, 4ug; Detection:Catalog No:109878 1:1000) with mouse brain tissue lysate 2640ug.