

# GST Tag antibody



Catalog Number: 111178

## Product name

GST Tag antibody

## Specificity

Tag

## Antibody description

GST Tag Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein.

## Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: FN315687). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

## Validations



Western blot of GST-tagged fusion protein with anti-GST-tag (Catalog No:117323) at various dilutions.