

LCTL antibody

Catalog Number: 112178

Product name

LCTL antibody

Specificity

Human; other species not tested.

Antibody description

LCTL Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 39 kDa, 65 kDa, 70 kDa

Preparation

This antibody was obtained by immunization of LCTL recombinant protein (Accession Number: NM_207338). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112178(LCTL Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:112178(LCTL Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG