

NHLRC1 antibody

Catalog Number: 113178

Product name

NHLRC1 antibody

Specificity

human; other species not tested.

Antibody description

NHLRC1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells. Positive IHC detected in human heart. Observed molecular weight by Western-blot: 45-47kd

Preparation

This antibody was obtained by immunization of NHLRC1 recombinant protein (Accession Number: NM_198586). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:113178(NHLRC1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human heart using Catalog No:113178(NHLRC1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:113178(NHLRC1 antibody) at dilution of 1:50 (under 40x lens)