

CRYBB1 antibody



Catalog Number: 107178

Product name

CRYBB1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CRYBB1 Mouse Monoclonal antibody. Positive WB detected in C6 cells, mouse eye tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of CRYBB1 recombinant protein (Accession Number: NM_001887). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



C6 cells were subjected to SDS PAGE followed by western blot with Catalog No:107178(CRYBB1 Antibody) at dilution of 1:1000