

# Anti-OR6K2 antibody



Catalog Number: 141578

## Product name

Anti-OR6K2 antibody

## Specificity

Human OR6K2

## Antibody description

Rabbit Polyclonal to Human OR6K2

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human OR6K2, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 0.5-5 µg/ml

## Validations



Items	Lanes	A	B
Sample (whole cell lysate)		HL60	HepG2
Sample Volume (µg/lane)		30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	0.5-5 µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (1:1), at 1:5000 dilution.		
Developed using Odyssey imaging system.			
Explanation	Predicted band size : 37 kDa Observed band size : 37 kDa		

Human OR6K2 Western blot (WB) 15197