

# Mecr antibody

Catalog Number: 112578

## Product name

Mecr antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

Mecr Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Positive IHC detected in human ovary tumor. Observed molecular weight by Western-blot: 65kd

## Preparation

This antibody was obtained by immunization of Mecr recombinant protein (Accession Number: NM\_025297). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112578(Mecr antibody) at dilution of 1:3000



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:112578(Mecr antibody) at dilution of 1:100 (under 10x lens)