

Anti-Adenosine A1 Receptor Rabbit antibody



Catalog Number: 168578

Product name

Anti-Adenosine A1 Receptor Rabbit antibody

Specificity

Human

Antibody description

Adenosine A1 Receptor Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Adenosine A1 Receptor

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

1:1000


1:50

1:20

1:20-1:100

1:20

Validations



Western blot detection of Adenosine A1 Receptor in K562,C6,3T3,HeLa cell lysates using Adenosine A1 Receptor Rabbit mAb(1:1000 diluted).Predicted band size:37kDa.Observed band size:37kDa.