

# BRCA2 antibody



Catalog Number: 108578

## Product name

BRCA2 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

BRCA2 Rabbit Polyclonal antibody. Positive IHC detected in rat brain tissue. Positive WB detected in human brain tissue, HeLa cells. Observed molecular weight by Western-blot: 384kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000059). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117225(BRCA2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded rat brain using Catalog No:117225(BRCA2 antibody) at dilution of 1:50 (under 10x lens)