

Human ANXA13/ANX13 recombinant protein



Catalog Number: 519578

General Information

Gene Name Synonym

Annexin A13, Annexin XIII, Annexin-13, Intestine-Specific Annexin, ISA, ANXA13, ANX13

Protein Construction

A DNA sequence encoding the Human ANXA13/ANX13 Gly2-His316 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 100mM NaCl, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ANXA13/ANX13 recombinant protein