

# MYO16 antibody



Catalog Number: 112978

## Product name

MYO16 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

MYO16 Rabbit Polyclonal antibody. Positive IHC detected in human bladder. Positive WB detected in human brain tissue, A431 cells. Observed molecular weight by Western-blot: 123kd,210kd, 106kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_015011). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112978(MYO16 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human bladder using Catalog No:112978(MYO16 antibody) at dilution of 1:50 (under 10x lens)