

# Anti-FGFRL1 antibody



Catalog Number: 175978

## Product name

Anti-FGFRL1 antibody

## Specificity

Human, Mouse

## Antibody description

Rabbit polyclonal antibody to FGFRL1

## Preparation

This antigen of this antibody was synthesized peptide derived from fgfr-5 . at aa range: 130-210

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1.; Western Blot analysis of 3T3 cells using FGFR-5 Polyclonal Antibody