

Anti-CD8 alpha Rabbit antibody



Catalog Number: 166978

Product name

Anti-CD8 alpha Rabbit antibody

Specificity

Human

Antibody description

CD8 alpha Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human CD8

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000


IHC: 1/50-1/200

ICC/IF: 1/50

Validations




Immunohistochemistry of CD8 in paraffin-embedded Human tonsil using CD8 Rabbit pAb at



dilution 1/100

Western blot detection of CD8 in THP-1 cell lysates using CD8 Rabbit pAb(1:1000 diluted).Predicted band size:26KDa.Observed



band size:29KDa.

Western blot detection of CD8 in K562,CHO-K1,HeLa cell lysates using CD8 Rabbit pAb(1:1000 diluted).Predicted band size:26KDa.Observed band size:29KDa.