

ELK1 antibody

Catalog Number: 110278

Product name

ELK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ELK1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in Jurkat cells, HepG2 cells, mouse pancreas tissue. Observed molecular weight by Western-blot: 62 kDa

Preparation

This antibody was obtained by immunization of ELK1 recombinant protein (Accession Number: XM_017029339). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110278(ELK1 antibody) at dilution of 1:800



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:110278(ELK1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)