

Anti-PWWP2B antibody



Catalog Number: 176278

Product name

Polyclonal

Anti-PWWP2B antibody

Ig Type

Specificity

Rabbit IgG

Human

Applications

Antibody description

WB, IHC-P

Rabbit polyclonal antibody to PWWP2B

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human pwwp2b. aa range:268-317

IHC-P: 1:100-1:300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Validations

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1:: Western Blot analysis of LOVO cells using PWWP2B Polyclonal Antibody