

# Anti-Cyclin H (1A7) Mouse antibody



Catalog Number: 168278

## Product name

Anti-Cyclin H (1A7) Mouse antibody

## Ig Type

Mouse IgG2b

## Specificity

Human, Mouse

## Applications

WB

## Antibody description

Cyclin H (1A7) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human Cyclin H protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of Cyclin H in Jurkat, K562, HeLa and 3T3 cell lysates using Cyclin H Mouse mAb (1:1000 diluted). Predicted band size: 38KDa. Observed band size: 38KDa.