

Human WDYHV1/C8orf32/NTAQ1 recombinant protein



Catalog Number: 519278

General Information

Gene Name Synonym

Protein N-terminal glutamine amidohydrolase,WDYHV1,Protein NH2-terminal glutamine deamidase,N-terminal Gln amidase,Nt(Q)-amidase,C8orf32, NTAQ1, REF: C1009

Protein Construction

A DNA sequence encoding the Human WDYHV1/C8orf32/NTAQ1 Met1-Cys205 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of PBS,100mM GSH,1% TritonX-100,15% glycerol,pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human WDYHV1/C8orf32/NTAQ1 recombinant protein