

# STARD13 antibody



Catalog Number: 115678

## Product name

STARD13 antibody

## Specificity

Human; other species not tested.

## Antibody description

STARD13 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HeLa cells. Observed molecular weight by Western-blot: 125 kDa

## Preparation

This antibody was obtained by immunization of STARD13 recombinant protein (Accession Number: NM\_001243466). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115678(STARD13 Antibody) at dilution of 1:1000