

NPY2R Polyclonal Antibody



Catalog Number: 163378

Product name

NPY2R Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to NPY2R

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human NPY2R (NP_000901.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - NPY2R Polyclonal Antibody

Western blot analysis of extracts of U-87MG cells, using NPY2R antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.