

# NUP50 antibody



Catalog Number: 113378

## Product name

NUP50 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NUP50 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 50kd,65kd

## Preparation

This antibody was obtained by immunization of NUP50 recombinant protein (Accession Number: XM\_017028537). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113378(NUP50 antibody) at dilution of 1:1000