

# TTC30A antibody



Catalog Number: 116378

## Product name

TTC30A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TTC30A Rabbit Polyclonal antibody. Positive WB detected in rat testis tissue, A431 cells, HEK-293 cells, mouse testis tissue. Observed molecular weight by Western-blot: 70-80 kDa

## Preparation

This antibody was obtained by immunization of TTC30A recombinant protein (Accession Number: NM\_152275). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



rat testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116378(TTC30A Antibody) at dilution of 1:1000