

Synaptotagmin-9 antibody



Catalog Number: 115778

Product name

Synaptotagmin-9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Synaptotagmin-9 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human kidney tissue, mouse eye tissue. Observed molecular weight by Western-blot: 52kd

Preparation

This antibody was obtained by immunization of Synaptotagmin-9 recombinant protein (Accession Number: BC029605). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115778(SYT9 antibody) at dilution of 1:500