

C6orf141 antibody



Catalog Number: 108778

Product name

C6orf141 antibody

Specificity

Human; other species not tested.

Antibody description

C6orf141 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of C6orf141 recombinant protein (Accession Number: BC036917). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108778(C6orf141 Antibody) at dilution of 1:1000