

# Human PCSK9/NARC1/PSEC0052 (AVI tag) recombinant protein



Catalog Number: 519778

## General Information

### Gene Name Synonym

Proprotein Convertase Subtilisin/Kexin Type 9, Neural Apoptosis-Regulated Convertase 1, NARC-1, Proprotein Convertase 9, PC9, Subtilisin/Kexin-Like Protease PC9, PCSK9, NARC1

### Protein Construction

A DNA sequence encoding the Human PCSK9/NARC1/PSEC0052 Gln31-Gln692 is expressed with a AVI tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 50mM HEPES, 150mM NaCl, pH 7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human PCSK9/NARC1/PSEC0052 (AVI tag)  
recombinant protein