

Human CMPK1 (His Tag) recombinant protein



Catalog Number: 502088

General Information

Gene Name Synonym

Deoxycytidylate kinase; Nucleoside-diphosphate kinase; Uridine monophosphate/cytidine monophosphate kinase

Protein Construction

A DNA sequence encoding the human CMPK1 (P30085) (Met1-Gly196) was fused with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Kinase activity untested

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human CMPK1 consists of 214 amino acids and has a calculated molecular mass of 24.5 kDa. The recombinant protein migrates as an approximately 27 kDa band in SDS-PAGE under reducing conditions.

Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, pH 7.4, 10% glycerol

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CMPK1 Protein (His Tag) SDS-PAGE