

Anti-Phospho-BTK (Tyr223) Rabbit antibody



Catalog Number: 167088

Product name

Anti-Phospho-BTK (Tyr223) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IP

Antibody description

Phospho-BTK (Tyr223) Rabbit polyclonal antibody

Dilutions

WB: 1/1000-1/2000

Preparation

Antigen: A synthesized peptide derived from human Phospho-BTK (Y223)

IP: 1/20

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Validations

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal



Western blot analysis of Phospho-BTK (Y223) expression in Raji cell lysate treated with pervanadate.