

Anti-Myc-Tag (9E10) Mouse mAb (HRP Conjugate)



Catalog Number: 169088

Product name

Anti-Myc-Tag (9E10) Mouse mAb (HRP Conjugate)

Ig Type

Mouse IgG1

Specificity

Tag

Applications

WB, ELISA, IHC-P, IP

Antibody description

Myc-Tag (9E10) Mouse mAb (HRP Conjugate)

Dilutions

ELISA: 1/100 - 1/500

Preparation


Antigen: Synthetic peptide corresponding to Human c-Myc aa 400 to the C-terminus (C terminal) conjugated to Keyhole Limpet Haemocyanin (KLH).

dilution.WB:1/10000 - 1/50000

Formulation

In PBS(pH 7.4) containing with 1mg/ml BSA 0.03% Proclin300 and 50% glycerol.

Validations



Storage

Ship at 4°C. Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot analysis of extracts from CHO-K1 cells, untransfected (A) or transfected with different FLAG-fusion proteins (B, C), using Myc-Tag (9E10) Mouse mAb (HRP Conjugate) Mouse mAb (1:1000 diluted).

Clonality

Monoclonal