

ETV4 antibody

Catalog Number: 110488

Product name

ETV4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ETV4 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 61 kDa

Preparation

This antibody was obtained by immunization of ETV4 recombinant protein (Accession Number: BC007242). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110488(ETV4 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:110488(ETV4 antibody) at dilution of 1:50 (under 40x lens)