

RB1 antibody



Catalog Number: 114488

Product name

RB1 antibody

Specificity

Human; other species not tested.

Antibody description

RB1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in Jurkat cells, A431 cells. Observed molecular weight by Western-blot: 106 kDa

Preparation

This antibody was obtained by immunization of RB1 recombinant protein (Accession Number: NM_000321). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:

WB: 1:200-1:1000

IF: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114488(RB1 Antibody) at dilution of 1:300

Immunofluorescent analysis of HepG2 cells using Catalog No:114488(RB1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)