

PIF1 antibody

Catalog Number: 113888

Product name

PIF1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PIF1 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 63-73 kDa

Preparation

This antibody was obtained by immunization of PIF1 recombinant protein (Accession Number: BC018978). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:113888(PIF1 antibody) at dilution of 1:800