

Anti-OR4D1 antibody



Catalog Number: 176188

Product name

Anti-OR4D1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR4D1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or4d1. aa range:261-310

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations



Fig1::; Western Blot analysis of various cells using Olfactory receptor 4D1 Polyclonal Antibody diluted at 1: 500