

# FAM40B antibody



Catalog Number: 110588

## Product name

FAM40B antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FAM40B Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, HeLa cells, HepG2 cells, Jurkat cells, K-562 cells, MCF-7 cells, mouse colon tissue, mouse testis tissue, Raji cells. Observed molecular weight by Western-blot: 95 kDa

## Preparation

This antibody was obtained by immunization of FAM40B recombinant protein (Accession Number: NM\_001134336). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:110588(FAM40B Antibody) at dilution of 1:600