

# Human Frizzled-10/FZD10 (CD350) (His Tag) recombinant protein



Catalog Number: 502588

## General Information

### Protein Construction

A DNA sequence encoding the human FZD10 (NP\_009128.1) (Met1-Gly161) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Ile 21

### Molecular Mass

The recombinant human FZD10 consists of 152 amino acids and predicts a molecular mass of 17.5 KDa. It migrates as an approximately 25.4 KDa

band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human Frizzled-10 / FZD10 Protein (His Tag)  
SDS-PAGE