

MYO9B antibody

Catalog Number: 112988

Product name

MYO9B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYO9B Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, HL-60 cells, Jurkat cells. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 230-250kd

Preparation

This antibody was obtained by immunization of MYO9B recombinant protein (Accession Number: BC018108). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:112988(MYO9B Antibody) at dilution of 1:600



IP Result of anti-MYO9B (IP:Catalog No:112988, 5ug; Detection:Catalog No:112988 1:500) with HepG2 cells lysate 2800ug.