

PNPLA3 Polyclonal Antibody



Catalog Number: 164988

Product name

PNPLA3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PNPLA3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 377-481 of human PNPLA3 (NP_079501.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - PNPLA3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PNPLA3 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.