

TAF6L antibody



Catalog Number: 115988

Product name

TAF6L antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TAF6L Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse thymus tissue, Raji cells. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of TAF6L recombinant protein (Accession Number: BC008785). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115988(TAF6L antibody) at dilution of 1:500