

DLEC1 antibody

Catalog Number: 109988

Product name

DLEC1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity


Human, Mouse; other species not tested.

IF: 1:20-1:200

Antibody description

DLEC1 Rabbit Polyclonal antibody. Positive IP detected in PC-3 cells. Positive WB detected in human testis tissue, A549 cells, PC-3 cells. Positive IF detected in PC-3 cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 80-90 kDa, 46 kDa

Validations




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_007335). Purification method: Antigen affinity purified.

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109988(DLEC1 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:

WB: 1:200-1:1000


Immunohistochemical of paraffin-embedded human testis using Catalog No:109988(DLEC1 antibody) at dilution of 1:50 (under 10x lens)



DLEC1 antibody


Catalog Number: 109988

Immunohistochemical of paraffin-embedded human testis using Catalog No:109988(DLEC1 antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-DLEC1 (IP:Catalog No:109988, 4ug; Detection:Catalog No:109988 1:300) with

PC-3 cells lysate 1100ug.



Immunofluorescent analysis of PC-3 cells using Catalog No:109988(DLEC1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)