

Anti-ALK-3 / BMPR1A antibody



Catalog Number: 100288

Product name

Anti-ALK-3 / BMPR1A antibody

Immunogen

[Human ALK-3 / BMPR1A \(His Tag\) recombinant protein](#)

Specificity

Human ALK-3 / BMPR1A

Antibody description

Rabbit polyclonal to ALK-3 / BMPR1A

Preparation

Produced in rabbits immunized with purified, recombinant Human ALK-3 / BMPR1A (rh ALK-3 / BMPR1A; NP_004320.2; Met1-Arg152). ALK-3 / BMPR1A specific IgG was purified by Human ALK-3 / BMPR1A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA:0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human ALK-3 / BMPR1A. The detection limit for Human ALK-3 / BMPR1A is < 0.039 ng/well.

IHC-P: 0.05-1 µg/ml

Validations



ALK-3 / BMPR1A Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry